

## Pulse Oximeter

Date: \_\_\_\_\_

Estimate Number: \_\_\_\_\_

Payment Terms: \_\_\_\_\_



### Product Features

- Name:** Bluetooth Pulse Oximeter
- Brand:** BalaMeng
- Packaging:** 1pcs/box, 100boxes/carton
- Intended Use:** Focus on blood oxygen testing for all ages. Blood oxygen pulse rate easy to detect with finger.
- Delivery:**
  - Gross weight: 5.45kg/carton.
  - Air volume weight: 4.85kg/carton.
  - Courier volume weight: 5.81kg/ carton.



Display screen:	OLED/LCD/ Digital display
SpO2 measuring range:	70~100%
SpO2 measurement accuracy:	±2% (70-100%) lower than 70% undefined
PR measuring range:	25~250BPM
PR measurement accuracy:	±2% or 2BPM
Resolution ratio:	168*80PX
Power dissipation:	lower than 30mA
Battery specifications:	Two AAA batteries
Size:	58*32*34mm



Quantity	Unit Price	Cost
	<b>Subtotal:</b>	
	<b>Tax:</b>	
	<b>Total:</b>	

