

Pulse Oximeter

Date: 2020.6.25

Estimate Number: _____

Payment Terms: _____



Display screen:	OLED/LCD/ Digital display
SpO2 measuring range:	70~100%
SpO2 measurement accuracy:	±2% (70-100%) lower than 70% undefined
PR measuring range:	25~250BPM
PR measurement accuracy:	±2% or 2BPM
Resolution ratio:	168*80PX
Power dissipation:	lower than 30mA
Battery specifications:	Two AAA batteries
Size:	58*32*34mm

Product Info:

- **Name:** Bluetooth Pulse Oximeter
- **Brand:** BalaMeng
- **Packaging:** 1pcs/box, 100boxes/carton
- **Intended Use:** Focus on blood oxygen testing for all ages. Blood oxygen pulse rate easy to detect with finger.
- **Delivery:** Gross weight 5.45kg/carton, air volume weight: 4.85kg/carton, courier volume weight: 5.81kg/carton

Features:

- **Small, lightweight, one-click measurement**

Simple and stylish appearance, rounded feel, detect blood oxygen anytime, anywhere

- **4 Seconds fast detection, industry leader**

Quick identification of physiological signal, only need 4 seconds.

- **Health data App**

The app displays blood oxygen (SpO₂), heart rate (PR), perfusion index (PI) in real time of the monitor, and save the measured data, which can be viewed and shared at any time.

- **Soft silicone finger membrane, comfortable to wear**


Special cavity structure allows silicone finger pad to adapt to the shape of finger.

- **Automatic shutdown in 10 seconds, continuous testing ≥24 hours**

Intelligent system identification and judgement, shut down if there is no measurement reaction for 10 seconds, intelligent low-power chip, continuous detection ≥ 24 hours

- OLED/LCD/Digital display

High-definition display with OLED/LCD/Digital display organic light-emitting secondary control parts. Can be displayed in 4 directions, multi-directional display.



Oximeter Certificates:

The following are the certificates of FDA, CE, RoHS:

- FDA certificate



FISCAL YEAR 2020 CERTIFICATE OF REGISTRATION

This Certifies that:

SHENZHEN MANRIDY TECHNOLOGY CO.,LTD

3rd Floor, Mingzesheng Building, Industrial Park Road, Dalang Street, Longhua District, Shen zhen, China

Was registered with US Food & Drug Administration, Center for Devices and Radio logical Health, pursuit to the Code of Federal Regulations 21 CFR 807. Such registration has been verified, with the registrant' s authorization, by TXB Rhineland Testing Services Corp Limited

Owner/Operator Number :10074234

Listing Number	Product Code(s)	Device Name	Proprietary Name
D403819	DQA	Oximeter	pulse oximeter SP001,SP002 ,MP010, MP020

TXB will confirm that such registration remains effective upon request and presentation of this certificate until the end of the calendar year stated above,unless said registration Is terminate after Issuance of this certification makes no other representations or warranties.nor does certificate make any representations or warrants to any person pr entity other than the named certificate holder.for whose sole benefit it ts issued.This certificate does not denote endorsement or approval of the certificate-holder' s device or establishment by the U.S. Food and Drug Administration.TXB assumes no affiliated to any person or entity In connection with We foregoing.

Pursuant to 21 CFR 807.39, " Registration of a device Establishment or assignment of a restoration number does not in any way denote approval.of the Establishment or its products ,Any reinterpretation mat create an Pinkerton of official approval because of registration or possession of a registration number is misleading and constitutes misleading" The US Food and Drug Administration does not Issue a certificate of registration nor does the Food and Drug Administration recognize a certificate of registration.TXB Is not affiliated with the U\$ Food and Drug Administration.



TXB Rhineland Testing Services Corp Limited

Web:HTTP://WWW.FDA.gov/

TXB Rhineland Testing Services Corp Limited

Issued:Apr.2,2020

Expiration date:Dec.31,2020

Cert.No:TXB20200519FDA05

E-mail: Lacy@xb-lab.com

WWW.xb-lab.com

● CE certificate

شهادة – 증명서 – Certificate – 證明書 – Сертификат

Shenzhen ZKT Technology Co., Ltd.
2F, Building 1, Maozhoushan Industrial Park, Houting Community,
Shajing Street, Bao'an District, Shenzhen, China



Certificate of Compliance

Certificate Number: ZKT-20200407979C

Certificate's Holder : SHENZHEN MANRIDY TECHNOLOGY CO.,LTD.
3rd Floor , Mingze Sheng Building, Industrial Park Road, Dalang Street, Longhua District, Shenzhen,China

Manufacturer : SHENZHEN MANRIDY TECHNOLOGY CO.,LTD.
3rd Floor , Mingze Sheng Building, Industrial Park Road, Dalang Street, Longhua District, Shenzhen,China

Trade Mark : **manridy**

Product : pulse oximeter

Model(s) : SP001, MP010, SP002, MP020

Test Standard : EN 55014-1:2017
EN 61000-3-2:2014, EN 61000-3-3:2013
EN 55014-2:2015

This Attestation of Compliance is issued on a voluntary basis for electrical equipment below the voltage limits of EMC Directive 2014/30/EU. The essential requirements are fulfilled accordingly based on the technical specifications applicable at the time of issuance. See also notes overleaf. It is only valid in connection with the test report number: ZKT-20200407979E.



This Certificate of Conformity is based on single evaluation of the submitted sample(s) of the above mentioned product. It does not imply an assessment of the whole product and relevant Directives to be observed.

● RoHS certificate

 **BST** | A RELIABLE TESTING FOR TRUST
GLOBAL TESTING AND CERTIFICATION PRECISION SERVICE CLOUD FACTORY

Certificate of Compliance

Certificate Number: BSTSH200514079801CC

Applicant : SHENZHEN MANRIDY TECHNOLOGY CO.,LTD.
3rd Floor , Mingze Sheng Building, Industrial Park Road, Dalang Street,
Longhua District, Shenzhen,China

Manufacturer : SHENZHEN MANRIDY TECHNOLOGY CO.,LTD.
3rd Floor , Mingze Sheng Building, Industrial Park Road, Dalang Street,
Longhua District, Shenzhen,China

Product Name : PULSE OXIMETER

Trade Name : **manridy**

Main Test Model : SP001

Additional Model : MP010, SP002, MP020

Test Standard : IEC 62321-4:2013+AMD1:2017, IEC 62321-5:2013,
IEC 62321-6:2015, IEC 62321-7-1:2015,
IEC 62321-7-2:2017, IEC 62321-8:2017

As shown in the Test Report No. : BSTSH200514079801CR

The EUT described above has been tested by us and found in compliance with the council RoHS 2 Directive 2011/65/EU Annex II (EU) 2015/863 as last amended by Directive (EU) 2017/2102. This certificate is only valid in conjunction with the test report.



RoHS


Tony Qian
Approved Signatory
May 19, 2020

Dongguan BST Testing Co., Ltd
Add: A1201-1204 Xinsanqi of Dongbao Road, Dongcheng District, Dongguan, Guangdong, China
Certificate Search: <http://www.bst-lab.com>, Tel: 400-8829628, 800-9990305, E-mail: christina@bst-lab.com